ITEM NO.:58608

3 Pagen Versand und Handelsges. mbH

St.-Jöris Strasse 16-28-52477 Alsdorf-Hoengen

GB

Including:

- 1* Foot Mat
- 1* Main Device
- 1* USB Charging Cable
- 1* User Manual

Suggest use one time per day.

Automatic turn off after 15 Seconds, if it is not used.

Power button is the plus button and to turn it of you must use the minus button.

- Use the plus or minus buttons to raise or lower the intensity of the 9 ranges
- Use the M-button to change the 6 different modes: 1. MASSAGE, 2. KNEAD, 3. SCRAPE, 4. ELBOW PRESSING, 5. ACUPUNCTURE, 6. CUPPING JAR
- Using time are 15 minutes
- Only one use per day
- Lowering the strength can prolong the battery life. The use time of battery power will be shortened by adjusting the strength.

People with these conditions CAN NOT use the foot mat:

- Anyone fitted with a pace maker or automatic implanted cardiac defibrillator(AICD)
- Anyone that is pregnant
- Anyone with deep vein thrombosis (DVT) to avoid mobilising the clot
- Anyone that have high blood pressure and anyone that have heart disease.

Please **DO NOT** use this product for the following parts, otherwise may cause an accident or skin, body issues.

- Chest, head (for example face)
- Wounds
- Recently operated body parts
- Metal and plastic implant site
- Spine
- tattoos

Do Not apply to body parts other than specified in this manual.

General Safety Information

This EMS Unit should never be used to mask or relieve undiagnosed pain.

Before treating any symptoms of pain, You should consult your physician.

For your safe and correct use of the device, it is recommended you to read the following guides before using it:

Precaution

- *Use a low intensity and shorter duration to introduce the muscles and avoid stimulation.
- *Keep out of reach of children and pets.
- *This product is not waterproof. Do not put it into the water or wash with water
- *Do not use the product if there is a sprain or wound in the feet, sole, ankle, achilles and calf.

Other important safety precautions.

1. Only use this device for the purpose of physical therapy,
2. Only use the accessories supplied with the foot mat.
3. Do not use this device in places with high humidity such as the bathroom.
4. Do not overload the electrical outlet.
5. If the device malfunctions, disconnect it from the power source and contact customer services, Do not attempt to repair the device by yourself. As do so will invalidate the warranty.
6. Do not let the device come in contact with any metal object, such as belt buckle or jewelry.

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The Enterprise Department Ltd.