Instruction Manual

When using this or any other electrical appliance, always follow these basic safety precautions:

⚠ DANGER

- Do not sit on pad.
- Do not use while sleeping.
- Never position pad underneath the body.
- Do not fold pad.
- Avoid sharp folds that may damage the product.

⚠ DANGER

Hot Surfaces.

DANGER: to reduce the risk of burns, electric shock and fire, this product must be used in accordance with the following instructions.

1. READ ALL INSTRUCTIONS.
2. DO NOT USE WHILE SLEEPING.
3. NEVER LEAVE THIS APPLIANCE UNATTENDED, ESPECIALLY IF CHILDREN ARE PRESENT.
4. DO NOT USE ON AN INFANT.
5. DO NOT USE ON AN ANIMAL.
6. THIS PAD IS NOT TO BE USED ON OR BY AN INVALID, SLEEPING OR UNCONSCIOUS PERSON, OR A PERSON WITH POOR BLOOD CIRCULATION UNLESS CAREFULLY ATTENDED.
7. THIS PAD IS NOT TO BE USED ON A PERSON WITH DIABETES.
8. DO NOT USE ON AREAS OF INSENSITIVE SKIN.
9. BURNS CAN OCCUR REGARDLESS OF CONTROL SETTING, CHECK SKIN UNDER PAD FREQUENTLY.
10. DO NOT USE IN AN OXYGEN ENRICHED ENVIRONMENT.
11. DO NOT USE PINS OR OTHER METALLIC MEANS TO FASTEN THIS PAD IN PLACE.
12. DO NOT SIT ON OR CRUSH PAD-AVOID SHARP FOLDS.
13. POSITION PAD ON TOP OF, AND NOT UNDERNEATH, THE BODY. HEAT LEVEL MAY INCREASE WHEN HEAT IS TRAPPED.
14. NEVER USE PAD WITH LINIMENTS, OINTMENTS OR SALVES, ESPECIALLY THOSE CONTAINING HEAT-PRODUCING INGREDIENTS, AS SKIN BURNS COULD RESULT.
15. NEVER PULL THIS PAD BY THE SUPPLY CORD.
16. NEVER UNPLUG PAD FROM ELECTRICAL OUTLET WITH A WET HAND.
17. IF USING AN EXTENSION CORD WITH THIS PAD, THE MARKED ELECTRICAL RATING OF THE CORD SET/EXTENSION CORD SHOULD BE AS GREAT AS ELECTRICAL RATING OF THE PAD (SEE ELECTRICAL SPECIFICATIONS SECTION). USE THIS PAD ONLY ON A 120 VOLT AC ELECTRICAL OUTLET. UNPLUG WHEN NOT IN USE.
18. DO NOT USE THE CORD AS A HANDLE.
19. CAREFULLY EXAMINE THE INNER VINYL PAD BEFORE EACH USE. DISCARD THE PRODUCT IF THE INNER VINYL PAD SHOWS ANY SIGN OF DETERIORATION.
20. DO NOT BEND OR PINCH THE CORD AND DO NOT USE THIS PAD IF IT HAS A DAMAGED POWER CORD OR PLUG, IF IT MALFUNCTIONS, HAS BEEN DROPPEN, OR IS DAMAGED. THIS PAD HAS NO USER-SERVICEABLE PARTS. IF FOR ANY REASON THIS PAD DOES NOT FUNCTION SATISFACTORILY, SEE WARRANTY FOR SERVICE INFORMATION.
21. DO NOT TAMPER WITH THIS PAD IN ANY WAY.

22. LOOP CORD LOOSELY WHEN STORING,TIGHT WARAPPING MAY DAMAGE CORD AND INTERNAL PARTS.
23. THIS PRODUCT IS INTENDED SOLELY FOR NON-COMMERCIAL,NON-INDUSTRIAL,PERSONAL USE. DO NOT USE THIS PRODUCT OUTDOORS OR FOR ANY OTHER PURPOSE.

The Safety Plug

This heating pad has a polarized plug (one blade is wider than the other) as a safety feature.

This plug will fit into a polarized outlet only one way. If the plug doesnot fit fully into the outlet, reverse the plug.

If it still does not fit, contact a qualified electrician.

DO NOT ATTEMPT TO DEFEAT THIS SAFETY FEATURE.

Care and Maintenance

1. Disconnect the power cord from the electrical outlet.

2. Make sure the Pad Connector is disconnected from the Control Connector by pressing the tabs on the sides of the Control Connector while pulling the connectors apart.

3. Hand wash the Pad using a mild detergent, cool water.

4. Line dry or lay the Pad flat to dry. Do not tumble dry.

5. Make sure the product is completely dry before the next use.

6. To store the product, loosely roll the Heating Pad toward the Storage Straps, loop the Straps over the rolled ends, and then store it in a cool, dry, dark location after each use. Avoid sharply folding the Heating Pad.

Electrical Specifications

Voltage	Frequency	Wattage	Amperage
AC220-240V	50Hz	70W	0.7A

Frequently Asked Questions

When should a heating pad be used?

Whenever a hot form of temporary relief is desired for personal comfort.

What temperature does this product reach?

Maximum temperature does not exceed the Underwriters Laboratories temperature limit of 176 °F (80°C) when set to "Hi". Temperature varies by heat setting;you may need to use a lower heat setting depending on the specific application and your personal preference.

How much time does it take to warm up to each temperature selection?

Temperature Selection	Indicator	Heat Up Elapsed Time
Low	1	6 min
Medium	2	8 min
High	3	10 min

DISPOSAL

- 3PAGEN uses reusable or recycled packaging wherever possible.

- Please dispose of all packaging, paper and cardboard in accordance with the recycling regulations of your local authority.

- At the end of the product's life, please enquire about disposal at an authorised household waste recycling centre in your municipality.

Safety instructions:

- Read the instructions for use and all safety instructions carefully before using the product.

- Use the product only for its intended purpose and in accordance with the manufacturer's instructions.

- Keep the product away from children: Many products are not suitable for children and should be kept out of the reach of children.

- Check the product regularly for damage or wear that could affect safety.

- Pay attention to all warnings and safety labelling on the product or packaging.

- Keep the product in a safe place to avoid accidents.

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